
BISCO

TheraBase™

Self-Adhesive Calcium Releasing Base/Liner

Instructions for Use

BISCO, Inc.
1100 W. Irving Park Road
Schaumburg, IL 60193
U.S.A.
1-847-534-6000
1-800-247-3368

IN-234R3
Rev. 1/21

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed Dentist.

en	TheraBase™*
Self-Adhesive Calcium & Fluoride Releasing Base/Liner	

GENERAL INFORMATION

TheraBase is a dual-cured, calcium and fluoride releasing, self-adhesive base/liner. **TheraBase** chemically bonds to tooth structure, and releases calcium and fluoride. It is radiopaque, allowing for easy identification on radiographs.

Intended user: Licensed dentist

Indications for Use:

- Base or Liner under restorations.

Warnings:

- Avoid contact with the skin; unpolymerized (meth)acrylate resins may cause skin sensitization in susceptible persons. In case of contact, wash skin with soap and water.
- Saliva and oral fluid contamination will seriously compromise dentin bonding.

Cautions:

Cross-contamination:

- Product may contain items that are designed for one-time use. Dispose of used or contaminated tips. Do not clean, disinfect or reuse.
- Use of commonly available hygienic protective covering to avoid any contamination of syringes during treatment is recommended. If syringe becomes contaminated, discard. Do not clean or disinfect.

Precautions:

- Working and setting times will be accelerated in the warm intra oral environment.
- Do not use hydrogen peroxide, EDTA, or eugenol containing materials prior to the adhesive procedure as these materials may inhibit bonding.
- Unpolymerized resins may cause skin sensitization in susceptible persons. Avoid contact with the skin. In case of contact with the skin, wash thoroughly with soap and water.
- To prevent polymerization or clogging, the Dual-Syringe static mixing tip should be left in place until the next application. If used intraorally, it should be wiped with disinfectant prior to storage.
- See individual component labels for specific expiration dates.
- Safety data sheet available on request.
- Safety data sheet available at www.bisco.com.

TECHNICAL INFORMATION

Working Time: Minimum 1 minute (includes mixing) at 22°C (71.6°F)

Setting Time: Maximum 5 minutes at 37°C (98.6°F)

Working time and setting time may vary based on storage conditions, temperature, humidity, etc.

INSTRUCTIONS FOR USE

Preparation Procedures

1. If a pulp exposure occurs, BISCO recommends using TheraCal LC as a pulp capping agent.
2. Clean the preparation with pumice and water. Rinse thoroughly.
3. Remove all surface water using a strong stream of air on the preparation for 3 to 5 seconds. Do not desiccate. (Protect the preparation from contamination. When applicable, use of a rubber dam is recommended.)

Mixing Procedures

1. Remove the syringe cap from the dual syringe.
2. Dispense a small amount of material onto a mixing pad to eliminate any voids in each chamber of the dual syringe.
3. Attach the mixing tip to the dual syringe by lining up the key and keyway. Then turn the locking mechanism clockwise to lock in place.
4. Press the plunger to mix and dispense the TheraBase. (See "TECHNICAL INFORMATION" section for working and setting time information.)

Basing/Lining Procedures

1. Apply TheraBase to the dentin surfaces of the prepared cavity directly from the dispensing syringe.
2. Light cure for 20 seconds or allow to self-cure. (When using self-cured mode, light curing the margins even briefly will expedite the final set.)

HELPFUL HINTS

- TheraBase left on a mixing pad may not appear to set. This is due to oxygen inhibition of the polymerization process. Rest assured the product will polymerize (set) in a normal fashion in the restoration.

HYGIENE: Use of commonly available hygienic protective covering to avoid any contamination of the TheraBase dual syringes during the treatment is recommended.

DISPOSAL: Refer to community provisions relating to waste. In their absence, refer to national or regional provisions relating to waste.

STORAGE: Store at room temperature (20°C/68°F - 25°C/77°F)

WARRANTY: BISCO, Inc. recognizes its responsibility to replace products if proven to be defective. BISCO, Inc. does not accept liability for any damages or loss, either direct or consequential, stemming from the use of or inability to use the products as described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

TO REORDER ACCESSORIES: X-81270P – 20-gauge auto-mix cannula tips (30 count).

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